


pEU-E01-MCS-TEV-GST-C1, C2 vectors

Multiple cloning site information



pEU-E01-MCS-TEV-GST-C1 (4358 bp)

EcoRV ↓ ↓ *XhoI* ↓ *BamHI* ↓ *SmaI* ↓ *KpnI* ↓ *NotI* ↓ *SpeI*
 GA TAT CTC GAG GAT CCC GGG TAC CGC GGC CGC ACT AGT
 Leu Glu Asp Pro Gly Tyr Arg Gly Arg Thr Ser

GAA AAC CTG TAT TTT CAG GGC	ATG GAA	GST	TCG GAT TAA
Glu Asn Leu Tyr Phe Gln Gly	Met Glu		Ser Asp STOP

↑ TEV cleavage site

pEU-E01-MCS-TEV-GST-C2 (4359 bp)

EcoRV ↓ ↓ *XhoI* ↓ *BamHI* ↓ *SmaI* ↓ *KpnI* ↓ *NotI* ↓ *SpeI*
 GAT ATC TCG AGG ATC CCG GGT ACC GCG GCC GCA CTA GTC
 Ile Ser Arg Ile Pro Gly Thr Ala Ala Ala Leu Val

GAA AAC CTG TAT TTT CAG GGC	ATG GAA	GST	TCG GAT TAA
Glu Asn Leu Tyr Phe Gln Gly	Met Glu		Ser Asp STOP

↑ TEV cleavage site